

NU BAND+AGE

THE MOST ADVANCED SECUREMENT DEVICE
PROVEN TO PROTECT, SUPPORT AND ASSIST
IN THE ACCELERATION OF WOUND HEALING.

{ NU BANDAGE TECHNOLOGY IS CHANGING THE FACE OF WOUND CARE. }



Easy to apply. Hand washable and reusable for same patient. Nu Bandage is available vented (2x10V-above) and vented padded (2x10VP-see back).

■ NO Trauma ■ NO Skin Tears ■ NO Waste

Nu Bandage is non-absorbent, easy to apply and reusable. It holds dressings, sub-bandages and medical devices of all types in place without the use of sticky tapes, messy adhesives, or clumsy fasteners. All while providing a protective second layer for those with skin sensitivities to pressure or touch.

IDEAL FOR TREATING:

- Pressure & venous stasis ulcers
- Diabetic wounds & ulcers
- Aids in the prevention of skin tears
- Patients who are sensitive to pressure & touch
- Excellent for burns & skin grafts
- Vascular & sensitive skin issues, especially in senior patients

PATIENT BENEFITS INCLUDE:

- Decreased tears vs. conventional bandages
- Easy access for improved medicine application
- Excellent airflow
- Reduced bruising of thin, easily traumatized skin
- Soft & gentle feel with a cooling effect
- Consistent, even pressure over wounds
- Fully articulating & contouring for greater patient comfort

NU BANDAGE IS:

- Available in vented & vented-padded options
- Easy to apply & wear
- Non-woven, non-absorbent, lightweight & soft
- Safe, biodegradable & latex-free
- Naturally anti-microbial
- Hand washable & reusable for same patient... NO buckles, fasteners or tapes


NU BAND+AGE

An Innovative Product from New Bandage Technology, LLC
WOUND CARE FOR THE 21st CENTURY

MADE IN THE USA
Patent Pending




Available in vented padded (2x10VP - above, available in white only) and vented (2x10V see front).
Easy to apply. Hand washable and reusable for same patient. No buckles, fasteners or tape necessary.



SUPERIOR IN RESULTS AND COST

Patients having suffered at least one skin tear in the previous 6 months who wore Nu Bandage experienced nearly **40%** fewer skin tears versus conventional alternatives. With a total cost savings of over **87%**.



See nubandage.com for references and details

{ Wasted spending in healthcare currently amounts to about \$765 billion, including \$105 billion in excessive prices and \$55 billion in missed PREVENTION opportunities. }

Product and Ordering Information:

Model (2" Series)	Inner Diameter x Length	Circumference Range	Recommended Use	Color
NB-2x5V (Vented)	5cm (2") ID x 12.5 cm (5") Length	20 - 40 cm	Arm/Foot/Hand/Lower Leg	White
NB-2x10V (Vented)	5cm (2") ID x 25.5 cm (10") Length	20 - 40 cm	Arm/Foot/Hand/Lower Leg	White
NB-2x10VP (Vented Padded)	5cm (2") ID x 25.5 cm (10") Length	20 - 35 cm	Arm/Foot/Hand/Lower Leg	White
NB-2x5S (Solid)	5cm (2") ID x 12.5 cm (5") Length	20 - 40 cm	Arm/Foot/Hand/Lower Leg	White
NB-2x10S (Solid)	5cm (2") ID x 25.5 cm (10") Length	20 - 40 cm	Arm/Foot/Hand/Lower Leg	White

Visit us at nubandage.com for more information and instructional video.

New Bandage Technology LLC
777 South Flagler Drive
Suite 800, West Tower
West Palm Beach, FL 33401



Scan QR Code with your Smartphone and watch our video!

Please contact your local distributor below for ordering and product information.

MEDIQUE MEDICAL SUPPLIES (MGBR) Inc.
5900, Andover Avenue
Town of Mount Royal, Quebec H4T 1H5
Tel: 514-342-4294 or 1-800-793-6210
Fax: 514-342-2294 or 1-800-793-4671
www.mediquemed.com

NU BAND+DAGE

An Innovative Product from New Bandage Technology, LLC
WOUND CARE FOR THE 21st CENTURY

MADE IN THE USA
Patent Pending