




PRESSURECARE

Pressure Relief & Patient Positioning



MEDIQUE
www.mediquemed.com
514-342-4294
1(800)-793-6210



Trulife Overview

Trulife is an internationally based and managed group engaged in the creation, development, manufacture and marketing of niche healthcare products. In addition to Pressurecare pads for the Operating Tables, the Group activities encompass product ranges, Orthopaedics, Prosthetics, Breastcare and Wheelchair cushions. For further details on all our product ranges please visit our website: www.trulife.com



Trulife products are distributed by more than 300 distributor companies worldwide. Our primary markets are the USA, Canada, the UK and Ireland, where we have our own dedicated marketing organisations.


The Group is headquartered in Dublin, Ireland. European manufacturing, research & development and marketing operations are conducted through wholly owned subsidiaries in Ireland and England. Orthopaedic clinical services are provided throughout England, Scotland and Wales.

In the Americas, Trulife trades through wholly owned subsidiary companies. Customer service and distribution facilities are based in Jackson, Michigan, Poulsville, Washington and Trenton, Ontario. Manufacturing and research & development of Orthopaedics, Prosthetics and Custom Breastcare products is conducted in Poulsville, Washington, Jackson, Michigan and Trenton, Ontario.

Mission Statement

The core philosophy is a belief in serving customers' needs through excellence in product innovation, quality in everything we do, extensive education and training offerings, and development of patient care through committed and dedicated people.

INTRODUCTION	2
Pressurecare	2
Trulife Ranges	3
PRODUCT GUIDE	4
Table Pad Protection	4
Head & Neck Protection	7
Arm & Elbow Protection	10
Body Protection	12
Leg & Knee Protection	17
Ankle, Heel & Foot Protection	18
Speciality Pad Protection	20
Paediatric Products	21
Accessories	25
POSITIONING GUIDE	26
Supine	26
Prone	27
Lateral	28
Lithotomy	29
TECHNICAL INFORMATION	30
Interface Pressure Maps	30
General Use	32
Index	33



[1] "In the UK Pressure Sores absorb around 4% of the National Healthcare expenditure"

[2] "2.5 million patients are treated for pressure sores in the USA"

Pressurecare Prevention of Peri-Operative Pressure Injury

Definition: Pressure Ulcer

'A pressure ulcer is localized injury to the skin and/or underlying tissue usually over a bony prominence, as a result of pressure, or pressure in combination with shear'
NPUAP / EPUAP 2011 [3]

Definition: Pressure Neuropathy

'In pressure neuropathy there is excessive pressure on a nerve root or a peripheral nerve leading to ischemia (not enough oxygen due to poor blood supply).. If the pressure continues, fibrosis tends to develop.' [4]

Who is at Risk?

All Operating Room (OR) patients due to level of Impaired Mobility.
Surgical patients still account for 20% to 30% of all pressure injuries.[5]
General Statistics on Prevalence of PUs
Germany (21.1%) [6], Australia (28.2%) [7], Italy (27%) [8]
US (11.5%) [9], Canada (16.5%) [9].

Areas Most at Risk

Bony Prominences – See Diagram
Head : Elbows : Sacrum : Heels : Shoulders : Hips : Buttocks : Knees.

Aetiology: What happens?

Prolonged Pressure ► Pressure on Blood Capillaries / Nerves
► Ischaemia & lack of Oxygen ► Skin Breakdown or Nerve Compression
► Pressure Ulcer / Neuropathy

Principles of Safe Positioning

▼ Level of prolonged pressure on bony prominences and nerve endings
► Avoid Circulatory and Respiratory Insufficiency
► Attain Optimal Surgical Position Safely
► Facilitate Role of Anaesthetist

Intervention - Pressure Injuries can be prevented! Use Preventative measures! Educate Staff!

Trulife offers an Extensive Range of Equipment for Pressurecare Management in the OR

The prevention of pressure injury assumes a fundamental role in health services economics and management. [10]



[1] Preventing Pressure Ulcers: The Goal Is Zero K. Duncan, R.N. The Joint Commission Journal on Quality and Patient Safety 2007; Vol 33 (10) 605 – 610

[2] The Cost of pressure ulcers in the UK G. Bennett, C. Dealey, J. Posnett age and ageing: Vol 33 (3) 230 – 235

[3] www.epuap.org www.npuap.org [4] www.medterms.com

[5] Patient Positioning: Education is key. www.surgicalproductsmag.com. 2010. Claire Williams Tager

[6] Lahmann N. et al., J. Ostomy/Wound Management.52: 20-33, 2006. [7] Gardner A. et al., J. Wound Practice and Research. 17:134-145, 2009.

[7] Lahmann N. et al., J. Ostomy/Wound Management.52: 20-33, 2006. [8] Capon A. et al., J. Advanced Nursing. 58: 263-272, 2007.

[9] Healthcare Management Council Inc. HMC: Hospitals lose \$2 Million in errors yearly. March 5, 2010. Assessed March 11,

2010 from: <http://www.hmc-benchmarks.com/2010/03/hmchospitals-lose-2-million-in-errors-yearly/>

[10] ECONOMIC EVALUATION OF PREVENTION OF PRESSURE ULCERS IN INTENSIVE CARE UNITS

Andrade, Melanie 1*, Nogueira, Fernanada 2, Morais, Ernesto 3: EPUAP Annual Meeting 2011

The incidence of OR-acquired pressure sores can be significantly reduced by using the most effective pressure relieving and positioning devices available and by employing a few basic concepts in positioning care. Education of OR staff is paramount.

The Trulife Pressurecare Ranges were developed as a result of our commitment to provide improved levels of pressure care products to both patients and healthcare systems. Vast amounts of in-depth research have gone into the development of each range. The Trulife Pressurecare range offers a simple and cost effective solution with a comprehensive offering to cover all surgical procedures.

Why select the Trulife Range?

- Extensive range
- Different budgets considered
- Medical Grade Silicone
- Latex Free
- MRI Compatible
- X-Ray Translucent
- Biocompatible
- Reusable
- Cost Saving
- Made in Ireland



Oasis

OASIS - THE ORIGINAL SILICONE SOLUTION

The Oasis range is an extensive range of reusable 100% silicone gel pads. The silicone we use is softer than both skin and underlying tissue. This means that pressure is relieved by the slight movement of the silicone gel and resulting dissipation of pressure/force across the product. It is the first range of products launched by Trulife in the OR field and remains a market leader in gel pads and positioners.

Oasis elite

OASIS ELITE - LIGHTWEIGHT & PRE CONTOURED

The Oasis Elite range combines silicone gel and light weight foam with a precontoured surface to provide superior pressure relief and patient protection in the OR. Anatomical pre-contouring facilitates greater conformity between gel and body shape; weight is then redistributed over a larger area. Oasis Elite products are more than 50% lighter than comparable gel-only products. They are designed for stability, easier handling and manipulation by medical personnel during preparation and surgery.

AZURE

AZURE - REDEFINING PRESSURECARE

The Azure range is the most advanced and innovative range of products available in the OR today. The Azure range offers the best pressure redistribution. Each product consists of two independent protective silicone layers that work in combination allowing the moulded products to effectively conform to the individual user's shape in a more contoured and less structured manner, resulting in effective protection against nerve damage and skin breakdown. The new approach by Trulife was to create a product that replicates and acts as an additional layer of protective tissue and skin. Azure has set a whole new standard of protection in the OR.

Surgical Positions products are suitable for:

Supine **S** Prone **P** Lithotomy **L** Lateral **LL** Paediatrics **P** * Indicates use

S P L LL



P

AZURE PROFILE PADS

The Azure Table Pads combine the ultra-soft and lightweight silicone material on the base and the cured silicone gel on top, with a unique squared profile that cushions the body during surgery. The pads are up to 20% lighter than comparative gel pads.

CODE	DESCRIPTION	DIMENSIONS
AZ610	Universal Table Pad	520 x 520 x 10mm
AZ611	¾ Length Table Pad	1150 x 520 x 10mm

S P L LL



P

AZURE TABLE PADS

The Azure Table pads can be used to protect the body in all surgical positions. The Azure Table Pads combine the ultra - soft and lightweight silicone material of the base and the cured silicone gel on top. The pads are up to 20% lighter than comparative gels.

CODE	DESCRIPTION	DIMENSIONS
AZ600	Operating Table Pad	520 x 520 x 10mm
AZ601	Operating Table Pad	1150 x 520 x 10mm

S P L LL



P

OASIS ELITE TABLE PADS

The Oasis Elite Pads can be used to protect the body in all surgical positions. The ultra light weight pad consists of a foam base and a gel top, allowing for easy carrying and manipulation by OR staff.

CODE	DESCRIPTION	DIMENSIONS
EL030	Universal Table Pad	520 x 520 x 30mm
EL031	¾ Length Table Pad	1150 x 520 x 20mm

S P L LL




P

OASIS STANDARD OPERATING TABLE PADS

These Oasis table pads provide protection over the entire surface area of the operating table and are recommended for use with patients in all positions. The pads are flexible, soft, comfortable and versatile and distribute weight evenly over the entire surface area.

CODE	DESCRIPTION	DIMENSIONS
OA030	Universal Table Pad	520 x 520 x 10mm
OA031	¾ Length Table Pad	1150 x 520 x 10mm
OA032	Full Length Table Pad	1800 x 520 x 10mm

S P L LL



P

OASIS PROFILE PAD

The Oasis Profile Pad has a series of 30mm² raised Trugel profiles that increase the overall surface area of the pad by 33%. This unique design allows for greater weight distribution and is particularly effective during extended procedures.

CODE	DESCRIPTION	DIMENSIONS
OA041	Universal Profile Pad	520 x 520 x 12mm
OA042	³ / ₄ Length Profile Pad	1145 x 520 x 12mm

S P L LL



P

OASIS 3 PIECE OPERATING TABLE PAD

The Oasis 3 Piece Operating Table Pad is designed for use on sectioned tables. The set consists of a head section, torso section and a foot section, with the latter two having a perineal cutout. The product provides a high degree of flexibility, versatility and protection during a variety of surgical procedures, particularly in the areas of Orthopaedics and Urology. Can also be used as a standard full length table pad.

CODE	DESCRIPTION	DIMENSIONS
OA036	Head Section	225 x 500 x 10mm
OA037	Torso Section	990 x 500 x 10mm
OA038	Foot Section	500 x 500 x 10mm
OA039	3 Piece Set	-



P

OASIS FULL LENGTH PAEDIATRIC PAD

Used to protect paediatric patients during all surgical procedures, this pad can also be used as a mattress overlay in hospital beds.

CODE	DESCRIPTION	DIMENSIONS
OA033	Paediatric Full Length Pad	605 x 335 x 10mm

S P L LL

**OASIS HIP PAD**

The Oasis Hip Pad is a flat pad specifically designed for use on hip fracture table tops.

CODE	DESCRIPTION	DIMENSIONS
OA035	Hip Pad	520 x 520 x 10mm

P

S P L LL

**OASIS PERINEAL TABLE PAD**

The Oasis Perineal Pad provides excellent protection to the sacral area, a point of maximum pressure. It is essential for all patients in the lithotomy position.

CODE	DESCRIPTION	DIMENSIONS
OA040	Perineal Table Pad	520 x 520 x 10mm

P

**INCUBATOR PAD**

This unique ultra soft foam pad provides a soft and comforting base for neonates while providing all-round protection from pressure sores and nerve damage.

CODE	DESCRIPTION	DIMENSIONS
OA300	Incubator Pad	620 x 340 x 12mm

P

S L LL



AZURE HEAD PROTECTOR

The dual layer Azure Head Protector is innovatively designed to contour to the individual patient's head, conforming comfortably while providing effective protection to the patient.

CODE	DESCRIPTION	DIMENSIONS
*AZ100	Infant/Adolescent/Adult	260 x 250 x 50mm



P * Infant/Adolescent only

P



OASIS ELITE PRONE HEAD REST

Developed in conjunction with the Operating Staff at a leading hospital, the Oasis Elite Prone Head Rest is anatomically designed so that the patient's head can comfortably rest face down with respiratory tubes in place. The unique design provides excellent support, stability, comfort and protection during surgery in the prone position, with openings on both sides. The product is lightweight, durable and easy to handle.

CODE	DESCRIPTION	DIMENSIONS
*EL142	Adolescent/ Adult	273 x 220 x 135mm
EL 141	Adult	280 x 240 x 140mm



P * Adolescent only

S



OASIS ELITE SUPINE HEAD REST

This Supine Head Rest is ideal for positioning the patient during eye surgery and other supine procedures. Ergonomically shaped to position the patient's head, while effectively protecting against the onset of pressure sores.

CODE	DESCRIPTION	DIMENSIONS
*EL152	Adolescent	275 x 230 x 65mm
EL151	Adult	275 x 237 x 90mm



P

* Adolescent use only

S

L

LL



OASIS ELITE CLOSED HEAD RING

Unique gel & foam combination conforms comfortably to the patient's head, this head ring is suitable for use in many surgical procedures.

CODE	DESCRIPTION	DIMENSIONS
*EL011	Infant/Adolescent	Ø140 x 34mm
EL012	Adult	Ø200 x 50mm




P

* Infant/Adolescent only

P

LL



OASIS ELITE OPEN HEAD RING

Combination of gel and foam provides stability and protection to the patient's face during surgery. An opening allows for intubation.


CODE	DESCRIPTION	DIMENSIONS
*EL021	Infant/Adolescent	Ø140 x 34mm
EL022	Adult	Ø200 x 54mm



P

* Infant/Adolescent only

S L LL



P

OASIS CLOSED HEAD RING

Available in a wide variety of size options, the Oasis Closed Head Ring protects the head, neck and ears during all types of surgery and is suitable for use in many procedures.

CODE	DESCRIPTION	DIMENSIONS
OA010	Neonate/New born	Ø90 x 20mm
OA011	Infant / Adolescent	Ø140 x 35mm
OA012	Adolescent/ Adult	Ø200 x 45mm



P LL



P

OASIS OPEN HEAD RING

Soft, comfortable and durable, the Oasis Open Head Ring supports the patient's face, while allowing easy access for intubation during surgery.

CODE	DESCRIPTION	DIMENSIONS
OA020	Neonate/New born	Ø90 x 22mm
OA021	Infant / Adolescent	Ø140 x 35mm
OA022	Adolescent/ Adult	Ø200 x 50mm



S L LL



P

OASIS FACE PAD

The Oasis Face Pad is ergonomically designed to protect the patient's forehead from pressure caused by endotracheal tubes.

CODE	DESCRIPTION	DIMENSIONS
OA140	Face Pad	165 x 80 x 8mm

S P LL



P


AZURE ARM PROTECTOR

The Azure Arm Protector assumes the shape of the individual patient's arm, acting as a second layer of tissue and offering excellent protection against nerve damage.

CODE	DESCRIPTION	DIMENSIONS
AZ200	Infant/Adolescent/Adult	330 x 140 x 30mm
AZ201	Adolescent/ Adult	460 x 180 x 35mm



S P L LL



OASIS ELITE CONTOURED ARM BOARD PAD

This anatomically shaped armboard pad is designed to provide added support and stability to the arm in a variety of procedures.

CODE	DESCRIPTION	DIMENSIONS
EL101	Adult	510 x 150 x 35mm
EL102	Adolescent/ Adult	395 X 140 X 40mm



S P L LL



P

OASIS ELITE ARMBOARD PAD

This flat armboard pad has a slight indentation to cradle the arm and provide effective pressure relief during various surgical procedures.

CODE	DESCRIPTION	DIMENSIONS
EL072	Adolescent/Adult	600 x 125 x 20mm
EL070	Infant/Adolescent/Adult	400 x 125 x 20mm

S P L LL




OASIS CONTOURED ARMBOARD PAD

This contoured armboard pad has been designed for use on a standard table extension. Its contoured shape provides excellent protection as well as keeping the arm in place.

CODE	DESCRIPTION	DIMENSIONS
OA101	Armboard Pad	510 x 150 x 35mm



S P L LL



OASIS ARMBOARD PADS

Oasis Armboard Pads provide the necessary protection whenever the arm is extended on an armboard. Suitable for all types of armboards.

CODE	DESCRIPTION	DIMENSIONS
OA070	Infant/Adolescent/Adult	395 x 115 x 10mm
OA071	Infant/Adolescent/Adult	600 x 200 x 10mm
OA072	Adolescent / Adult	600 x 125 x 10mm
OA073	Adolescent / Adult	605 x 427 x 10mm



P

S L LL




OASIS ULNAR / BRACHIAL NERVE PROTECTORS

The Oasis Ulnar / Brachial Nerve Protectors are essential for providing protection to the sensitive areas of the elbows and arms.

CODE	DESCRIPTION	DIMENSIONS
OA060	Adolescent/ Adult	380 x 115 x 8mm
OA061	Adolescent/ Adult	460 x 210 x 8mm
OA062	Adolescent/ Adult	460 x 380 x 8mm

P * Adolescent only

S



OASIS ARM RETAINER PAD

The Oasis Arm Retainer Pad has been specially designed to fit perspex / acrylic arm supports. It can be easily fitted with Velcro attachments.

CODE	DIMENSIONS
OA100	170 x 110 x 18mm
OA102	157 x 128 x 15mm



P

S L



AZURE SACRAL PROTECTOR

The Sacral area is an area of serious risk if left unprotected. The Azure Sacral protector has been designed to provide support in the supine and lithotomy positions by relieving pressure and focussing on peripheral support.

CODE	DESCRIPTION	DIMENSIONS
*AZ300	Adolescent/Adult	400 x 380 x 25mm
AZ301	Adult	450 x 445 x 30mm



P * Adolescent only

S P LL



AZURE UNIVERSAL POSITIONER

This versatile and lightweight product can be used to position and protect the patient in variety of surgical procedures.

CODE	DESCRIPTION	DIMENSIONS
AZ500	Adolescent/ Adult	425 x 155 x 60mm



P * Adolescent only

P LL



OASIS ELITE PRONE LATERAL POSITIONING SYSTEM

The Prone/Lateral Positioning System from Trulife is the result of extensive research into the needs of Operating Room personnel for patients during surgical procedures in the prone/lateral position.

Designed with both patient and surgical staff in mind, this system offers the ultimate pressure relieving protection and effective positioning all in one.



CODE	DESCRIPTION	DIMENSIONS
EL905	Oasis Elite Prone/Lateral Positioning System	675 x 500 x 100mm

REPLACEMENT PARTS AVAILABLE:

EL906	Upper Left Section /Lower Right Section	275 x 250 x 100mm
EL907	Upper Right Section /Lower Left Section	275 x 250 x 100mm
EL908	Mid Left Section /Mid Right Section	125 x 250 x 100mm
EL912	Base Plate	500 x 500 x 6mm

P



OASIS ELITE PRONE POSITIONER

The Prone Positioner gel and foam combination has been designed to provide effective positioning to the patient during surgical procedure in the prone position. For best protection, it is recommended to be used in conjunction with the Prone Head Rest.

CODE	DESCRIPTION	DIMENSIONS
EL903	Adult	610 x 530 x 130mm



P



OASIS ELITE LARGE PRONE POSITIONER

The Large Prone Positioner eliminates pressure on the abdomen and facilitates easy ventilation and circulation during long cases in prone position, i.e surgeries to the backs of legs and laminectomy procedure.

CODE	DESCRIPTION	DIMENSIONS
EL904	Adult	1115 x 520 x 140mm



LL



OASIS ELITE LATERAL POSITIONER

This positioning device permits a safe, quick and easy transition from supine to lateral. The patient's lower shoulder fits into the integral channel to provide exceptional patient comfort and nerve protection. The Lateral Positioner is designed for use in hip replacement, shoulder, kidney and other procedures. The built-in headrest and auxiliary support assists in proper spinal alignment and surgical site access.

CODE	DESCRIPTION	DIMENSIONS
*EL900	Adolescent/ Adult	720 x 500 x 160mm



P

* Adolescent only: Lateral Position

S P



OASIS ELITE DOME POSITIONER

This lightweight positioner can be used to position and protect the patient during various surgical procedures in all positions. It can be used in many different applications.

CODE	DESCRIPTION	DIMENSIONS
*EL210	Adolescent / Adult	350 x 150 x 75mm

P

* Adolescent only: For Prone Position

S P L LL



OASIS ELITE FLAT BOTTOM CHEST ROLLS

Chest Rolls provide excellent support and stability. The unique gel and foam combination means they are up to 75% lighter than comparative gel products, allowing for easy manipulation and handling by OR staff.

CODE	DESCRIPTION	DIMENSIONS
EL214	Adult	300 x 140 x 110
EL215	Adult	400 x 140 x 110
EL216	Adult	505 x 140 x 110
EL217	Adult	405 x 150 x 150
EL218	Adult	510 x 150 x 150

S P L LL



OASIS FLAT BOTTOM CHEST ROLLS

These 100% silicone gel chest rolls provide excellent positioning and stability and are available in a wide variety of sizes.

CODE	DESCRIPTION	DIMENSIONS
*OA208	Infant	150 x 80 x 50mm
*OA209	Infant	200 x 80 x 50mm
*OA210	Infant / Adolescent	245 x 80 x 50mm
OA211	Adult	305 x 95 x 75mm
OA212	Adult	405 x 95 x 75mm
OA213	Adult	505 x 92 x 75mm



P *Specific for Paediatrics

S P L LL



OASIS PATIENT POSITIONING BELT

The Oasis Patient Positioning Belt ensures patient safety by providing stability and tissue protection on the Operating Table. By using the clip-lock attachment the belt can be easily adjusted to fit all operating tables. In case of emergency, there is an easy release clip lock that can be removed quickly.

CODE	DIMENSIONS
OA110	505 x 85 x 8mm

S P L LL



P

OASIS EXTREMITY BELT

The Oasis Extremity Belt is used in many procedures for extra security in maintaining the position of the extremities on the operating table.

CODE

OA113

DIMENSIONS

315 x 45 x 10mm

S P L LL



P

OASIS UNIVERSAL PATIENT POSITIONERS

These 100% silicone positioners are suitable for safe and effective positioning during all surgical procedures.

CODE	DESCRIPTION	DIMENSIONS
OA050	Infant/Adolescent/Adult	150 x 60 x 40mm
OA051	Infant/Adolescent/Adult	330 x 70 x 50mm
OA052	Adolescent / Adult	290 x 100 x 60mm
OA053	Adolescent / Adult	345 x 120 x 70mm

S



OASIS PERINEAL POST PAD

The Oasis Perineal Post Pad reduces pressure on the patient in the groin area and is attached quickly and easily with Velcro. Its universal shape and size allows it to be used on most operating room fracture tables.

CODE

OA180

DIMENSIONS

300 x 260 x 10mm

S P LL



AZURE LEG PROTECTOR

The Azure Leg Protector provides overall protection to the heel and leg area. Slightly raised behind the knee, the leg protector offers protection to the popliteal nerve and conforms superbly around the leg. Overall positioning, protection and alignment are the key features of this unique product.

CODE	DESCRIPTION	DIMENSIONS
*AZ400	Adolescent / Adult	515 x 180 x 40mm



P * Adolescent only

L



AZURE CRUTCH STIRRUP PADS

These Crutch Stirrup Pads were designed for use with various types of legholders to protect the nerves and blood vessels in the leg during lithotomy procedures. Lightweight and flexible, the dual layer Azure pads come in pairs and protect against skin shear and nerve injury. Supplied in pairs.

CODE	DIMENSIONS
AZ410	370 x 205 x 25mm

LL



OASIS ELITE LATERAL LEG POSITIONER

This product is used to position the legs in a lateral position for kidney surgery and hip surgery. In addition to comfort, it facilitates optimal surgical positioning. This positioner also allows for knee and hip flexion. To use: Place patient on side, place upper limb on top (gel part of the product), and place lower limb underneath ledge.

CODE	DESCRIPTION	DIMENSIONS
EL800	Adult	650mm x 400mm x 230mm

P



OASIS ELITE KNEE PAD

The Knee pad is used for various procedures in the prone position to protect the knees against pressure damage and the onset of pressure sores.

CODE	DESCRIPTION	DIMENSIONS
EL401	Adult	510 x 150 x 35mm
*EL402	Adolescent/ Adult	400 x 150 x 35mm



P * Adolescent only

L



OASIS CRUTCH STIRRUP PADS

The Oasis Crutch Stirrup Pads are designed to fit most types of stirrups. Supplied in pairs, they come in a variety of sizes.

CODE	DESCRIPTION	DIMENSIONS
OA090	Adolescent / Adult	320 x 220 x 10mm
OA091	Adolescent / Adult	395 x 245 x 10mm
OA092	Adolescent / Adult	520 x 300 x 10mm



S

P



AZURE HEEL PROTECTOR

The heel area is one of extreme risk during any surgical procedure and special attention should be given to avoid any unnecessary damage. These lightweight protectors provide heel protection as well as stability and support. Supplied in pairs.

CODE	DESCRIPTION	DIMENSIONS
AZ401	Infant/Adolescent/ Adult	190 x 175 x 30mm



S

P



OASIS ELITE HEEL PADS

The Oasis Elite Heel Pads provide support and maximum protection for the heels, even during the longest surgical procedures. They support the Achilles and effectively remove pressure from the heels. Supplied in pairs, the products are lightweight and easy to use.

CODE	DESCRIPTION	DIMENSIONS
*EL130	Adolescent / Adult	180 x 100 x 70mm



* Adolescent only

S



P * Adolescent only

OASIS HEEL PADS

These 100% silicone gel pads are suitable for use in any supine procedure where the heels are at risk of skin breakdown and damage.

CODE	DESCRIPTION	DIMENSIONS
*OA130	Adolescent/Adult	180 x 100 x 70mm

S



P

OASIS HEEL CUPS

Oasis Heel Cups are designed to provide pressure relief for patient's heels. They are suitable for use during and after most surgical procedures and can be used also for patients who are confined to bed or are in a day chair.

CODE	DESCRIPTION	DIMENSIONS
OA131	Neonate/ Infant/ Adolescent/Adult	105 x 90 x 43mm



S L LL



P * Adolescent only: For Lateral

OASIS HEEL & ANKLE PROTECTOR

These unique pads are designed to fit around the heel and ankle like a boot with Velcro attachment. They provide protective padding for the heel when the patient is in the supine or lateral positions. Can also be used as a protector for the ankle and foot when using conventional ankle strap stirrups and lithotomy poles.

CODE	DESCRIPTION	DIMENSIONS
*OA200	Adolescent/ Adult	250 x 120 x 10mm



L



OASIS BOOT STIRRUP PAD

These pads provide excellent pressure relief and protection to the calf and heel and are designed for use on most types of boot stirrups. The pads are secured to the boot stirrup via Velcro.

CODE	DESCRIPTION	DIMENSIONS
*OA160	Adult	605 x 140 x 5mm

* Also suitable for Fowler's position

P



OASIS RELTON FRAME PADS

Designed to fit over the four support pads of the Relton frame and provide excellent protection to the patient when positioned for spinal surgery.

CODE	DESCRIPTION	DIMENSIONS
OA120	Adult	190 x 170 x 18mm




P



OASIS WILSON FRAME PADS

Wilson Frame Pads are designed to provide pressure relief for patients positioned on a Wilson Frame, used to maintain patients in a flexed spine position.

CODE	DESCRIPTION	DIMENSIONS
OA221	Adult	735 x 255 x 10mm




INCUBATOR PAD
OA300



AZURE PROFILE PADS
AZ610 / AZ611




AZURE TABLE PADS
AZ600 / AZ601




OASIS ELITE TABLE PADS
EL030 / EL031



OASIS STANDARD OPERATING
TABLE PADS
OA030 / OA031 / OA032




OASIS PROFILE PAD
OA041 / OA042



OASIS 3 PIECE OPERATING
TABLE PAD
OA036 / OA037 / OA038
/ OA039



OASIS FULL LENGTH
PAEDIATRIC PAD
OA033



OASIS PERINEAL TABLE PAD
OA040



AZURE HEAD PROTECTOR
AZ100



OASIS ELITE PRONE
HEAD REST
EL142



OASIS ELITE SUPINE
HEAD REST
EL152



OASIS ELITE CLOSED
HEAD RING
EL011



OASIS ELITE OPEN
HEAD RING
EL021



OASIS CLOSED
HEAD RING
OA010 / OA011



OASIS OPEN HEAD RING
OA020 / OA021



OASIS FACE PAD
OA140



AZURE ARM PROTECTOR
AZ200




OASIS ELITE ARMBORAD PAD
EL070



OASIS ARMBORAD PADS
OA070 / OA071 / OA072 / OA073



OASIS ULNAR / BRACHIAL
NERVE PROTECTORS
OA060 / OA061 / OA062



OASIS ARM RETAINER PAD
OA100 / OA102




AZURE SACRAL PROTECTOR
AZ300



AZURE UNIVERSAL POSITIONER
AZ500



OASIS ELITE LATERAL
POSITIONER
EL900




OASIS ELITE DOME
POSITIONER EL210



OASIS FLAT BOTTOM CHEST
ROLLS
OA208 / OA209 / OA210



OASIS EXTREMITY BELT
OA113



OASIS UNIVERSAL PATIENT
POSITIONERS
OA050 / OA051 / OA052
/ OA053



AZURE LEG PROTECTOR
AZ400




OASIS ELITE KNEE PAD
EL402



AZURE HEEL PROTECTOR
AZ401



OASIS ELITE HEEL PADS
EL130



OASIS HEEL PADS
OA130



OASIS HEEL CUPS
OA131



OASIS HEEL & ANKLE
PROTECTOR
OA200



OA005

Trulife
First Aid
Repair Kit

PRESSURECARE
Technical Guide

OEACARD12

Pressurecare for Operating Room, Decubitus and Nerve Injury Prophylaxis

- Original Gel Solution
- 100% Silicone
- Similar Texture to Human Skin

- Gel & Foam
- Pre-Contoured
- 50% Lighter than Oasis

- Anatomical Design
- Dual-Layer of Silicone
- 20% Lighter than Oasis

- Medical Grade Silicone
- Latex Free
- Reusable
- MRI Compatible
- Easy to Clean
- Biocompatible
- 3-ling Translucent
- Most Extensive Range

www.trulife.com

Distributor:

Technical Guide

Operating Room Poster

PROTECTION FOR ALL SURGICAL PROCEDURES

SUPINE POSITION

- Administration of General Anesthesia
- Abdominal Surgery
- Open Heart Surgery
- Surgery on Face, Neck or Mouth
- Most Surgery on Extremities

PRONE POSITION

- Surgery on the Back, Spine and the back of the legs

LATERAL POSITION

- Surgery of the Chest or Lungs
- Surgery on the Pelvis
- Hip Surgery

LITHOTOMY POSITION

- Urological Surgery
- Gynecological Surgery

www.trulife.com

Positioning Guide

SUPINE POSITION


INCLUDES TRENDLENBURG & REVERSE TRENDLENBURG

The patient is on their back with face toward ceiling, legs are uncrossed and arms are at the side or on arm boards.

- Administration of General Anesthesia
- Abdominal Surgery
- Open Heart Surgery
- Surgery on Face, Neck or Mouth
- Most Surgery on Extremities

Body Area	Product Description	Colors
HEAD & NECK	Head Protector	AZ000
	Shoulder Head Rest	EL000 / EL100 / EL200 / EL300 / EL400 / EL500
	Support Head Rest Pad	EL000 / EL100 / EL200 / EL300 / EL400 / EL500
ARM & ELBOW	Arm Restraints	AZ000 / AZ100
	Armboard Pads	EL000 / EL100 / EL200 / EL300 / EL400 / EL500
	Upper / Backrest Protector	AZ000 / AZ100
BODY	Torso Pads	AZ000 / AZ100
	Spinal Protector	EL000 / EL100 / EL200 / EL300 / EL400 / EL500
	Universal Protector	EL000 / EL100 / EL200 / EL300 / EL400 / EL500
LEG & KNEE	Leg Protector	AZ000
	Leg Pads	EL000 / EL100 / EL200 / EL300 / EL400 / EL500
	Foot & Ankle Protector	EL000 / EL100 / EL200 / EL300 / EL400 / EL500

www.trulife.com



BODY AREA	PRODUCT DESCRIPTION	CODES
HEAD & NECK	Head Protector Closed Head Rings Supine Head Rest Face Pad	AZ100 EL011 / 12 OA010 / 11 / 12 EL151 / 52 OA140
ARM & ELBOW	Arm Protectors Arm/Kneeboard Pads Ulnar / Brachial Protector Arm Retainer Pad	AZ200 / 201 EL070 / 72 / 101 / 102 OA101 / 070 / 071 / 072 / 073 OA060 / 061 / 062 OA100 / 102
BODY	Table Pads Sacral Protector Universal Positioner Dome Positioner Patient Positioners Positioning Belt Extremity Belt Perineal Post Pad Chest Rolls	AZ 600 / 601 / 610 / 611 EL030 / 31 OA030 / 31 / 32 / 35 / 36 / 37 / 38 / 39 / 40 / OA041 / 42 AZ300 / 301 AZ500 EL210 OA050 / 51 / 52 / 53 OA110 OA113 OA180 EL214 / 215 / 216 / 217 / 218 OA211 / 212 / 213
LEG & KNEE	Leg Protector	AZ400
ANKLE, HEEL & FOOT	Heel Protector Heel Pads Heel Cups Heel & Ankle Protector	AZ401 EL130 / OA130 OA131 OA200



PRONE



BODY AREA	PRODUCT DESCRIPTION	CODES
HEAD & NECK	Head Protector Open Head Rings	AZ100 EL021 / 22 OA020 / 21 / 22
	Prone Headrest	EL141 / 142
ARM & ELBOW	Arm Protectors Armboard Pads	AZ200 / 201 EL070 / 72 / 101 / 102 OA101 / 070 / 071 / 072 / 073
BODY	Table Pads	AZ 600 / 601 / 610 / 611 EL030 / 31 OA030 / 31 / 32 / 35 / 36 / 37 / 38 / 39 / 40 / OA041 / 42
	Prone Positioning System Prone Positioner Large Prone Positioner Universal Positioner Chest Rolls	EL905 / 906 / 907 / 908 / 912 EL903 EL904 AZ500 EL214 / 215 / 216 / 217 / 218 OA211 / 12 / 13
	Dome Positioner Patient Positioners Extremity Belt Positioning Belt	EL210 OA050 / 51 / 52 / 53 OA113 OA110
LEG & KNEE	Arm / Knee Pad Leg Protector	EL401 / 402 AZ400
ANKLE, HEEL & FOOT	Dome Positioner	EL210
SPECIALITY PAD	Relton Frame Pads	OA120
	Wilson Frame Pads	OA221



BODY AREA	PRODUCT DESCRIPTION	CODES
HEAD & NECK	Open Head Ring	EL021 / 022
	Head Protector	OA020 / 021 / 022
	Closed Head Rings	AZ100
	Face Pad	EL011 / 12 OA010 / 11 / 12 OA140
ARM & ELBOW	Arm/Knee Protectors	AZ200 / 201
	Armboard Pads	EL070 / 72 / 101 / 102
	Ulnar / Brachial Protector	OA101 / 070 / 071 / 072 / 073 OA060 / 061 / 062
BODY	Table Pads	AZ600 / 601 / 610 / 11 EL030 / 31 OA030 / 31 / 32 / 33 / 35 / 36 / 37 / 38 / 39 / 40 / OA041 / 42
	Lateral Positioner	EL900
	Prone Positioning System	EL905 / 906 / 907 / 908 / 912
	Universal Positioner	AZ500
	Patient Positioners	OA050 / 51 / 52 / 53
	Positioning Belt	OA110
	Extremity Belt	OA113
	Flat Bottom Chest Rolls	EL214 / 215 / 216 / 217 / 218 OA211 / 212 / 213
LEG	Azure Leg Protector	AZ400
ANKLE, HEEL & FOOT	Heel & Ankle Protector	OA200
	Oasis ELite Lateral Leg Positioner	EL800





BODY AREA	PRODUCT DESCRIPTION	CODES
HEAD & NECK	Head Protector Closed Head Rings	AZ100 EL011 / 12 OA010 / 11 / 12
	Face Pad	OA140
ARM & ELBOW	Arm Protectors Arm/Legboard Pads	AZ200 / 201 EL070 / 71 / 72 / 73 / 101 / 102 OA101 / 070 / 071 / 072 / 073
	Ulnar / Brachial Protector Arm Retainer Pad	OA060 / 061 / 062 OA100 / 102
BODY	Table Pads	AZ600 / 601 / 610 / 11 EL030 / 31 OA030 / 31 / 32 / 33 / 35 / 36 / 37 / 38 / 39 / 40 / OA041 / 42
	Incubator Pad Sacral Protector Patient Positioners Positioning Belt Extremity Belt Flat Bottom Chest Rolls	OA300 AZ300 / 301 OA050 / 51 / 52 / 53 OA110 OA113 EL214 / 215 / 216 / 217 / 218 OA208 / 209 / 210 / 211 / 212 / 213
LEG & KNEE	Crutch Stirrup Pad	AZ410 OA090 / 91 / 92
ANKLE, HEEL & FOOT	Heel & Ankle Protector Boot Stirrup Pad	OA200 OA160



WITHOUT PROTECTION





HEAD









ARMS



HEELS



Tested in accordance with protocol for interface pressure management Ref: FSAT.0017



“Any positioning device used to reduce pressure on patient's body parts and prevent tissue damage should have the documented ability to reduce capillary interface pressures to 32mmHg or less (5)”

⁵ Hendrick, J. & Thomson et al. "Pressure reduction products: making appropriate choices" Journal of Enterostomal Nursing 20 (Nov/Dec 1993)

Each and every product in the Trulife Pressurecare ranges have been subjected to the most rigorous of durability tests to ensure that they will withstand normal day to day handling in even the busiest surgical environment. In the design phase, tests are conducted in relation to the every day effects of compression, rolling, pulling and tearing.

The Trulife Pressurecare ranges consist of the most extensive range of gel pads and positioners available and can be used for a variety of applications in practically every surgical procedure.

The products are manufactured under strict ISO: 13485:2003 standards and carry the CE mark.

HOW CAN THE PRODUCTS BE CLEANED?

All products can be easily cleaned using standard OR environment detergents. It is not advised to soak the products over long periods. The products are not suitable for use in conjunction with an autoclave.

CAN THE PRODUCTS BE USED IN DIRECT CONTACT WITH THE PATIENT'S SKIN?

All products are latex-free and non-allergenic. The gel used in the products is a cured silicone and does not flow. It is an inert material and does not support the growth of bacteria.

ARE YOUR PADS X-RAY TRANSLUCENT?

Our pads are x-ray translucent, however, higher radiation levels may be needed.

CAN THE PRODUCTS BE HEATED OR COOLED?

The products can be cooled to -12°C without any adverse effects. As adhesion of velcro can be weakened by low temperatures, products with velcro should be allowed to return to room temperature prior to use.

The products can also be heated to a maximum of 40°C if required and can be safely used in conjunction with heated blankets.

CAN THE PADS BE USED WITH MRI EQUIPMENT?

The materials used in our Oasis/Elite/Fusion and Azure products are a combination of PU film and silicone. They contain no ferrous material and so are considered MRI friendly.

WHAT IF THE PRODUCT GETS CUT?

While a repair kit is available, in the event of the product being cut, it is highly recommended that the product be discarded to avoid the risk of infection.

DO THE PRODUCTS CARRY A GUARANTEE?

All products are guaranteed against manufacturing defects for a period of 2 years.

PLEASE NOTE:

- Due to our manufacturing process a slight variation in colour may occur between individual products.
- Product dependent a tolerance of +/- 10mm is applied to all product dimensions.
- The position of the model in certain images is intended to display the products only and may not be surgically accurate in all cases.
- We offer a guideline to help select products for different surgical positions. Please refer to good clinical practices when selecting products.

CODE	DESCRIPTION	DIMENSIONS (mm)	PAGE	CODE	DESCRIPTION	DIMENSIONS (mm)	PAGE
TABLE PAD PROTECTION				BODY PROTECTION			
AZ610	Azure Universal Table Pad	520 x 520 x 10	4	AZ300	Azure Sacral Protector Small	400 x 380 x 25	12
AZ611	Azure ¾ Length Table Pad	1150 x 520 x 10	4	AZ301	Azure Sacral Protector Large	450 x 445 x 30	12
AZ600	Operating Table Pad	520 x 520 x 10	4	AZ500	Azure Universal Positioner	425 x 155 x 60	12
AZ601	Operating Table Pad	1150 x 520 x 10	4	EL903	Oasis Elite Prone Positioner	610 x 530 x 130	13
EL030	Oasis Elite Universal Table Pad	520 x 520 x 30	4	EL904	Oasis Elite Large Prone Positioner	1115 x 520 x 140	14
EL031	Oasis Elite ¾ Length Table Pad	1150 x 520 x 20	4	EL905	Oasis Elite Prone Positioning System	675 x 500 x 100	13
OA030	Oasis Universal Table Pad	520 x 520 x 10	4	EL906	Upper Left Section/Lower Right Section	275 x 250 x 100	13
OA031	Oasis ¾ Length Table Pad	1150 x 520 x 10	4	EL907	Upper Right Section/Lower Left Section	275 x 250 x 100	13
OA032	Oasis Full Length Table Pad	1800 x 520 x 10	4	EL908	Mid Left Section/Mid Right Section	125 x 250 x 100	13
OA041	Oasis Universal Profile Pad	520 x 520 x 12	5	EL912	Base Plate	500 x 500 x 6	13
OA042	Oasis ¾ Length Profile Pad	1145 x 520 x 12	5	EL900	Oasis Elite Lateral Positioner	720 x 500 x 160	14
OA036	Oasis Head Section	225 x 500 x 10	5	EL210	Oasis Elite Dome Positioner	350 x 150 x 75	14
OA037	Oasis Torso Section	990 x 500 x 10	5	EL214	Oasis Elite Flat Bottom Chest Rolls	300 x 140 x 110	15
OA038	Oasis Foot Section	500 x 500 x 10	5	EL215	Oasis Elite Flat Bottom Chest Rolls	400 x 140 x 110	15
OA039	Oasis 3 Piece Set		5	EL216	Oasis Elite Flat Bottom Chest Rolls	505 x 140 x 110	15
OA033	Oasis Paediatric Full Length Pad	605 x 335 x 10	5	EL217	Oasis Elite Flat Bottom Chest Rolls	405 x 150 x 150	15
OA035	Oasis Hip Pad	520 x 520 x 10	6	EL218	Oasis Elite Flat Bottom Chest Rolls	510 x 150 x 150	15
OA040	Oasis Perineal Table Pads	520 x 520 x 10	6	OA208	Oasis Flat Bottom Chest Rolls	150 x 80 x 50	15
OA300	Incubator Pad	620 x 340 x 12	6	OA209	Oasis Flat Bottom Chest Rolls	200 x 80 x 50	15
HEAD & NECK PROTECTION				OA210	Oasis Flat Bottom Chest Rolls	245 x 80 x 50	15
AZ100	Azure Head Protector	260 x 250 x 50	7	OA211	Oasis Flat Bottom Chest Rolls	305 x 95 x 75	15
EL142	Oasis Elite Prone Head Rest Medium	273 x 220 x 135	7	OA212	Oasis Flat Bottom Chest Rolls	405 x 95 x 75	15
EL141	Oasis Elite Prone Head Rest Large	280 x 240 x 140	7	OA213	Oasis Flat Bottom Chest Rolls	505 x 92 x 75	15
EL152	Oasis Elite Supine Head Rest Adolescent	275x 230 x 65	8	OA110	Oasis Patient Positioning Belt	505 x 85 x 8	15
EL151	Oasis Elite Supine Head Rest Adult	275 x 237 x 90	8	OA113	Oasis Extremity Belt	315 x 45 x 10	16
EL011	Oasis Elite Closed Head Ring Adolescent	Ø140 x 34	8	OA050	Oasis Universal Patient Positioners	150 x 60 x 40	16
EL012	Oasis Elite Closed Head Ring Adult	Ø200 x 50	8	OA051	Oasis Universal Patient Positioners	330 x 70 x 50	16
EL021	Oasis Elite Open Head Ring Adolescent	Ø140 x 34	8	OA052	Oasis Universal Patient Positioners	290 x 100 x 60	16
EL022	Oasis Elite Open Head Ring Adult	Ø 200 x 54	8	OA053	Oasis Universal Patient Positioners	345 x 120 x 70	16
OA010	Oasis Closed Head Ring Paediatric	Ø90 x 20	9	OA180	Oasis Perineal Post Pad	300 x 260 x 10	16
OA011	Oasis Closed Head Ring Adolescent	Ø140 x 35	9	LEG & KNEE PROTECTION			
OA012	Oasis Closed Head Ring Adult	Ø 200 x 45	9	AZ400	Azure Leg Protector	515 x 180 x 40	17
OA020	Oasis Open Head Ring Paediatric	Ø90 x 22	9	AZ410	Azure Crutch Stirrup Pads	370 x 205 x 25	17
OA021	Oasis Open Head Ring Adolescent	Ø140 x 35	9	EL800	Oasis Elite Lateral Leg Positioner	650 x 400 x 230	17
OA022	Oasis Open Head Ring Adult	Ø200 x 50	9	EL401	Oasis Elite Knee Pad	510 x 150 x 35	17
OA140	Oasis Face Pad	165 x 80 x 8	9	EL402	Oasis Elite Knee Pad	400 x 150 x 35	17
ARM & ELBOW PROTECTION				OA090	Oasis Crutch Stirrup Pads	320 x 220 x 10	17
AZ200	Azure Arm Protector Short	330 x 140 x 30	10	OA091	Oasis Crutch Stirrup Pads	395 x 245 x 10	17
AZ201	Azure Arm Protector Long	460 x 180 x 35	10	OA092	Oasis Crutch Stirrup Pads	520 x 300 x 10	17
EL101	Oasis Elite Contoured Armboard Pad	510 x 150 x 35	10	ANKLE, HEEL & FOOT PROTECTION			
EL102	Oasis Elite Contoured Armboard Pad	395 x 140 x 40	10	AZ401	Azure Heel Protector	190 x 175 x 30	18
EL072	Oasis Elite Armboard Pad Long	600 x 125 x 20	10	EL130	Oasis Elite Heel Pads	180 x 100 x 70	18
EL070	Oasis Elite Armboard Pad Short	400 x 125 x 20	10	OA130	Oasis Heel Pads	180 x 100 x 70	19
OA101	Oasis Contoured Armboard Pad	510 x 150 x 35	11	OA131	Oasis Heel Cups	105 x 90 x 43	19
OA070	Oasis Armboard Pads	395 x 115 x 10	11	OA200	Oasis Heel & Ankle Protector	250 x 120 x 10	19
OA071	Oasis Armboard Pads	600 x 200 x 10	11	OA160	Oasis Boot Stirrup Pad	605 x 140 x 5	20
OA072	Oasis Armboard Pads	600 x 125 x 10	11	SPECIALTY PAD PROTECTION			
OA073	Oasis Armboard Pads	605 x 427 x 10	11	OA120	Oasis Relton Frame Pads	190 x 170 x 18	20
OA060	Oasis Ulnar/Brachial Nerve Protectors (S)	380 x 115 x 8	11	OA221	Oasis Wilson Frame Pads	735 x 255 x 10	20
OA061	Oasis Ulnar/ Brachial Nerve Protectors (L)	460 x 210 x 8	11	CONCLUSION			
OA062	Full Length	460 x 380 x 8	11	Paediatric Products			21
OA100	Oasis Arm Retainer Pad	170 x 110 x 18	11	Accessories			25
OA102	Oasis Arm Retainer Pad	157 x 128 x 15	11	Positioning Guide			26
				Interface Pressure Maps			30
				General Use			32

Advice to ensure long life and best Service from your Trulife Pressurecare products.

Thank you for choosing our Pressurecare product range.

Our products are used to aid the prevention of pressure ulcers and correct positioning in the surgical environment. A Patient's predisposition to this type of injury is based on many factors including intrinsic factors such as general health, age, weight and pre-existing medical conditions. The risk of patients acquiring pressure related injuries can be minimized by properly using the Trulife range of products, along with good clinical practices. Appropriate care and attention will ensure the long life and durability of your Trulife Pressurecare products.

General use of your Trulife Pressurecare products.

- In line with good clinical practice, assess the patient's condition and areas of potential risk prior to deciding which pads and positioners to use.
- Please ensure the patient is correctly positioned in a safe and secure manner at all times and ensure that no excess pressure is applied to any one area. Pay particular attention to areas of high risk.
- A sheet may be placed over the product; however, this will reduce its shear relieving properties. If a sheet is used, it should fit loosely over the products and be un-starched and free of wrinkles or folds.
- Products should be handled with care. Avoid placing sharp instruments on the products as they will puncture the film.
- Pads may be heated to body temperature in an oven, pre-heated to 40°C, or on a heated mattress. Please follow instructions of the oven or mattress manufacturer as appropriate. (Please refer to good clinical practices)
- The products can be cooled to -12°C without any adverse effects.
- Avoid pulling table pads while patients are lying on them, as excessive force will adversely affect them.
- Avoid pulling products along rough or sharp surfaces.
- If any damage occurs (e.g. cut, crack or tear) the product should be discarded. A repair kit is available but should only be used as a temporary solution.
- Do not expose the product to a naked flame.
- Use correct lifting techniques when lifting heavy products.

Oasis Extremity Belt (OA113) and Oasis Positioning Belt (OA110)

This product is designed to assist in stabilising the patient when positioned on the operating table, and is not designed to support weight of the patient.

Prone Positioners EL903 and EL904

- For Maximum protection, these products should be used in conjunction with the prone headrest (EL141), and a low chest roll or universal Positioner to elevate the feet.
- Ensure the patient is positioned safely and securely.
- Please ensure straps are securely wrapped around the side rails of the table.

Oasis Prone / Lateral Positioning System (EL905)

EL905: For Prone

- The EL905 is a versatile and flexible products that can be adjusted to suit the patient. For correct assembly please align black dots as indicated in Diagram 1 page 2
- For Maximum protection, the products should be used in conjunction with the prone headrest (EL141), the knee Pad (EL101/102) and a low chest roll, or universal Positioner to elevate the

- feet. Ensure the patient is positioned safely and securely. (Please refer to first image on page 3)
- To best maintain the adhesion of the Hook and Loop, peel the Positioner away from the base plate using the inside corner.
- Please ensure straps are securely wrapped around the side rails of the table.

EL905: For Lateral

- Set up the EL905 for Prone Position (See Image Diagram 1, page 2)
- Rotate the top two gel pieces, you will see one continuous black line (See Diagram 2, Page 2).
- Move the other 4 pieces of the EL905 down to the next Velcro piece, this will allow the EL101 to fit.
- The EL905 should be adjusted as per image, once in this position it should not be adjusted.
- Using the EL905 in Lateral, please use in conjunction with the EL022 / Open Head Ring and the EL101 Armboard Pad (Please see Diagram 3, page 2).
- An EL101 pad should also be used under the shoulder. (This will cover the black plastic base on the EL905 as illustrated in Diagram 2, page 2).The patient should never come in contact with the black plastic base of the EL905.
- Ensure the patient is positioned safely and securely.
- Please ensure straps are securely wrapped around the side rails of the table.

Advice to ensure long life and best Service from your Trulife Pressurecare products.

Cleaning of your Trulife Pressurecare products.

- Products should be disinfected before and after use.
- Products can be easily cleaned using standard operating room cleaning agents. Avoid soaking products in cleaning fluid.
- If product has Hook and Loop attachments keep it closed when not in use to reduce the build up of dust.

Handling

Gently remove pads from tables or surgical equipment. To transfer large table pads (not containing foam) you may roll the pad and carry it cradled in arms. To transport larger foam or gel products, fully support the pad with two hands or cradle it in arms when possible. Carts are recommended for transporting bulkier items.

Storage of your Trulife Pressurecare products

- Fold large table pads in sections of approximately 35cm.
- Avoid placing heavy loads or sharp objects on any of the products.
- Store products at room temperature.

Oasis Elite Table Pads (EL030 and EL031)

For best results store product with Gel side facing outwards.

Disposal instructions

Any disposal practice must be in compliance with all local and national laws and regulations. All packaging materials used for your Pressurecare product can be recycled.

In the unlikely event should any problems arise with your Pressurecare product, please contact your local distributor.

If returning a damaged product, please disinfect prior to sending.

If you have any questions regarding our product range or their applications, or have any comments regarding any Products, please contact your local distributor.

NOTE: PLEASE RETAIN THESE INSTRUCTIONS FOR USE, FOR REFERRAL.



MEDIQUE
www.mediquemed.com
514-342-4294
1(800)-793-6210