


# Samilast

## NET BANDAGES



This tubular synthetic and elastic net can secure bandages on any part of the body and that can replace traditional gauze bandages.



### BODY PARTS

- 0 Fingers and toes
- 1 Hand, wrist, foot, ankle
- 2 Arm, leg, elbow, knee
- 3 Arm, leg, elbow, knee, shoulder
- 4 Head, thigh, hip spica
- 5 Hip spica, thigh, chest, head
- 6 Perineum, axilla, chest, pelvis
- 7 Perineum, axilla, chest, pelvis
- 8 Perineum, axilla, chest, pelvis
- 10 Perineum, axilla, chest, pelvis

### ADVANTAGES :

- It can be cut in both the lengthwise and in the width direction without fraying which will hold dressings in difficult to dress and small areas.
- It can be used for fast and easy retention of all types of wound dressings.
- It allows easier access to inspect wounds which allows them to be handled quickly and efficiently.
- It is impervious against blood, sweat and even ointments.
- Most importantly, it allows the wound to be aired.

**COMPOSITES:** 38% polyester  
35% polyamide  
27% rubber

**MEDIQUE**, EXCLUSIF DISTRIBUTOR IN NORTH AMERICA.

### Technical data for Samilast products:

Size	Product code	Body parts	Stretched length	Unstretched width
0	SAM-S050M	Fingers and toes	50 m/54,7 vg	0,8 cm/0,3"
1	SAM-S150M	Hand, wrist, foot, ankle	50 m/54,7 vg	1,5 cm /0,6"
2	SAM-S250M	Small arm, leg, elbow, knee	50 m/54,7 vg	1,8 cm/0,7"
3	SAM-S350M	Large arm, leg, elbow, knee, shoulder	50 m/54,7 vg	2,0 cm/0,8"
4	SAM-S450M	Small head, thigh, hip spica	50 m/54,7 vg	2,3 cm/0,9"
5	SAM-S550M	Large head thigh, hip spica and very small chest	50 m/54,7 vg	2,6 cm/1,0"
6	SAM-S625M	Small perineum, axilla, chest, pelvis	25 m/27,3 vg	5,2 cm/2,0"
7	SAM-S725M	Medium perineum, axilla, chest, pelvis	25 m/27,3 vg	5,5 cm/2,2"
8	SAM-S825M	Large perineum, axilla, chest, pelvis	25 m/27,3 vg	6,3 cm/2,5"
10	SAM-S1025M	Extra large perineum, axilla, chest, pelvis	25 m/27,3 vg	10 cm/3,9"

Please note that data for stretched length was measured by a dynamometer with a 16kg tension.

The unstretched width can be slightly different due to surrounding factors such as heat, light, humidity and aging.

For more information about these products, please contact our customer service department:

# MEDIQUE



5900 Ave Andover  
Ville Mont-Royal, QC H4T 1H5



1-800-793-6210



1-800-793-4671



[customerservice@mediquemed.com](mailto:customerservice@mediquemed.com)



[www.mediquemed.com](http://www.mediquemed.com)